

Direct detection of AKT/PKB activity in a *Pten* knock out mouse model using dynamic peptide microarrays.

4046

ACR 2011



Riet Hilhorst¹, Liesbeth Houkes¹, Hanneke Korsten², Monique Mommersteeg¹, Jan Trapman² and Rob Ruijtenbeek¹

¹ PamGene International BV, Nieuwstraat 30, 5211 NL 's-Hertogenbosch, The Netherlands. Phone +31736158076, rihilhorst@pamgene.com, www.pamgene.com

² Department of Pathology, Josephine Nefkens Institute, Erasmus Medical Center, Rotterdam, The Netherlands.

Background

In many cancers, deregulation of kinases or phosphatases leads to activation of the AKT/PKB pathway. This pathway regulates many cellular processes including metabolism, differentiation, translation, proliferation and apoptosis, so is central in determining cell fate. In addition, this pathway is often hyperactive in tumor cells resistant to kinase inhibitor therapies (survival pathways). Several kinase inhibitors that influence AKT activity are investigated in the context of new pharmacotherapies in cancer. However, often a fraction of patients responds and it is not possible to identify those patients up front. It would be a major leap forward if response to therapy could be determined in an *in vitro* assay prior to the onset of therapy. Methods to directly measure this activity do not exist.

In this study, we explored:

- 1) Whether it is feasible to determine AKT activity profiles in lysates of hyperplastic prostates of PTEN knock-out mice (PSA-Cre;Pten-loxP/loxP mice).
- 2) Whether addition of AKT inhibitors affects the phosphorylation profiles in an inhibitor specific way.

PamChip® Cell-Based Kinase Assay




Fig. 1. A. PamChip® 96 microarray plate and peptide microarray used for kinase activity profiling. The solid support is three-dimensional and consists of pores in which the peptides are immobilized. B. Peptide phosphorylation through active kinases, or inhibition thereof, can be monitored in lysates of cell lines or tumor tissues by fluorescent readout, using a set of anti-phospho-ser/thr specific antibodies and a fluorescently labeled secondary antibody for detection. C. The signals per peptide are quantified and analysed using Bionavigator software.

Presence of Akt in *Pten* knockout mice




Fig. 2. The PSA-Cre;Pten-loxP/loxP mouse model. A. Effect of *Pten* knockout on histology of the prostate. B. Overexpression of pS473-AKT in hyperplastic cells of PSA-Cre;Pten-loxP/loxP mice. C. Presence of pS473-AKT, pT308-AKT, total AKT, MAPK and S6RP in prostate lysates. Equal amounts of protein were loaded on the gel.

Activity of Akt in *Pten* knockout mice




Fig 3. A. Comparison of phosphorylation profiles of lysates of PSA-Cre;Pten-loxP/loxP and wild type mice with recombinant AKT1. B. The data represented in a heatmap. Each horizontal line represents a peptide. Signals are the average of five replicates. C. Effect of inhibitors on signal intensity. The ratio of inhibited/non-inhibited is represented for AKT Inhibitors VIII and MK-2206 and PKC/PDK inhibitor UCN-01. Signals are the average of three replicates.

Western blots showed the presence of AKT in the wild type lysate and in the lysate of the hyperplastic prostate. Only the *Pten* knockout lysates contained AKT phosphorylated on S473 and T308, indicative of AKT activity (Fig. 2). Phosphorylation activity profiles revealed differential activity between PSA-Cre;Pten-loxP/loxP and wt mouse prostate lysates, with the *Pten* negative lysate having increased kinase activity in all biological replicates tested (Fig. 3A). Activity profiles differed from those of recombinant AKT1. In the *Pten* knockout lysates, addition of AKT inhibitors reduced signals on AKT substrate peptides and downstream peptides to the level observed in the lysates of the wt mice. Addition of non-AKT specific inhibitors like UCN-01 resulted in a different inhibition profile as compared to Akt inhibitors.

Methods

The PSA-Cre;Pten-loxP/loxP mouse prostate cancer model was described in [1] and [2]. Prostate tissues of three 4 months old PSA-Cre;Pten-loxP/loxP mice and three wt mice were lysed in MPER buffer with protease and phosphatase inhibitors. IHC and western blotting were performed using standard methods.

Lysates and recombinant full length AKT1 were incubated at least in triplicate on PamChip® peptide micro-arrays [3] comprising 256 ser/thr containing peptides derived from known human phosphorylation sites. Per array, 0.5 µg of protein was applied. As controls, incubations without ATP were performed. Peptide phosphorylation was monitored with a mixture of 6 anti-phosphoserine and anti-phosphothreonine antibodies. The effect of AKT inhibitors MK2206 and AKT Inhibitor VIII and non-AKT inhibitor UCN-01 (hydroxystaurosporine) on phosphorylation profiles was tested at 18 µM, 50 µM and 10 µM respectively.

References

1. Ma, X., Ziel-van der Made, A.C., Autar, B., van der Korput, H.A., Vermeij, M. et al (2005) Cancer Res 65 5730-5739
2. Korsten, H., Ziel-van der Made, A., Ma, X., van der Kwast, T. and Trapman, J. (2009) PLoS ONE 4(5): e5682.
3. Folkvord S et al. "Prediction of response to preoperative chemoradiotherapy in rectal cancer by multiplex kinase activity profiling.", Int J Radiat Oncol Biol Phys. 2010 Oct 1;78(2):555-62.

Conclusions

We measured AKT activity in a PTEN knock-out (PSA-Cre;Pten-loxP/loxP) mouse prostate tissue lysate, using a multiplex kinase assay. The activity of kinases in this tissue was increased compared to wild type tissue levels. In the presence of AKT and non-AKT inhibitors, differential modulation of kinase activity was observed.

The activity profiles of the PTEN/AKT pathway will support the studies on the complicated regulation of AKT activity. This new multiplex kinase activity assay will add a new dimension and a scope of applications like drug testing directly in patient derived tissue (biopsies). Detection of a-priori active resistance pathways in pre-dose biopsies is such an application with relevance to drug response prediction.