

# 4655

# In vitro tyrosine kinase activity profiling to identify molecular targets and predictive biomarkers in gastric cancer cell lines

AACR 2013

Yasuhiro Koh<sup>1</sup>, Masakuni Serizawa<sup>1</sup>, Rik de Wijn<sup>2</sup>, Riet Hilhorst<sup>2</sup>, Takako Nakajima<sup>3</sup>, Hirofumi Yasui<sup>4</sup>, Kazuyoshi Yanagihara<sup>5</sup>, Martijn Dankers<sup>2</sup>, Rob Ruijtenbeek<sup>2</sup>, Narikazu Boku<sup>3</sup>

<sup>1</sup>Drug Discovery and Development Division, Shizuoka Cancer Center Research Institute, Shizuoka, Japan, <sup>2</sup>PamGene International B.V., 's-Hertogenbosch, The Netherlands  
<sup>3</sup>Department of Clinical Oncology, St. Marianne University School of Medicine, Kawasaki, Japan, <sup>4</sup>Division of Gastrointestinal Oncology, Shizuoka Cancer Center, Shizuoka, Japan  
<sup>5</sup>Division of Genetics, National Cancer Center Research Institute, Japan



## Abstract

**Background:** Gastric cancer is among the leading causes of cancer-related deaths worldwide, notably in Asia. Though targeted therapies have been under clinical evaluation, the biological complexity of this cancer type may confer resistance to single-targeted therapies. Especially, scirrhous phenotype is known to be affected by multiple growth factors. To overcome this biological complexity, multi-targeted kinase inhibitors (MTKIs) have been clinically tested and were reported to be effective in a subset of patients. We wanted to elucidate the kinase activity signature of gastric cancer and to explore predictive biomarkers of MTKIs using gastric cancer cell lines by tyrosine kinase activity profiling.

**Methods:** The growth-inhibitory effect of the MTKIs sunitinib, sorafenib, pazopanib and the Akt inhibitor MK-2206 was evaluated by the MTT assay. 11 scirrhous and 14 non-scirrhous gastric cancer cell lines were grown until semi-confluent, lysed, aliquoted and stored at -80°C until use. The basal tyrosine kinase activity of the lysates was measured on PamChip® peptide micro-arrays, containing 144 peptides derived from known human phosphorylation sites. Quantification of peptide phosphorylation and data analysis were performed using BioNavigator and SIMCA P+ software.

**Results:** HSC-39, HSC-40A, KATO-III and HSC-43 scirrhous cells and SNU-16 non-scirrhous cells were sensitive to all three MTKIs with IC50 values < 1 μmol/L. Orthogonal partial least square-discriminant analysis (OPLS-DA) of basal tyrosine kinase activity profiles showed a clear segregation between sensitive cell lines and insensitive cell lines. Levels of putative target proteins such as PDGFR, VEGFR, Raf and Kit did not correlate with the sensitivity. All sensitive cell lines turned out to harbor an FGFR2 amplification suggesting that FGFR2 amplification may be a predictive biomarker of response to MTKI treatment. Kinases such as EGFR, EPHA, FGFR, IGF1R, INSR and Src contributed to this segregation suggesting that these are relevant to the sensitivity to MTKIs. MET overexpressing cell lines were resistant to treatment with the four inhibitors. Scirrhous cells were segregated from non-scirrhous cells by OPLS-DA among cell lines with signal-ring cell carcinoma cytology with R2, 0.994; Q2, 0.615. EGFR, ERBB2, 3, 4, EPHAs, FAK, IGF1R, INSR, MAP2Ks, RET and SRC seemed to be responsible for the segregation suggesting that they play an important role in scirrhous gastric cancer biology.

**Conclusions:** Kinase activity profiles reflect the biological complexity of gastric cancer. FGFR2 status correlates with sensitivity to sunitinib, pazopanib and sorafenib, whereas MET status does not. These data imply that tyrosine kinase activity profiling with PamChip® peptide arrays can be utilized to identify molecular targets and predictive biomarkers in gastric cancer and this should be further explored.

## Dynamic peptide micro-array technology




Figure 1. (A)(B) PamStation4 instrument and PamChip4 array plates used for kinase activity profiling. (C) The reaction takes place on spotted peptides on the array. Peptide phosphorylation through active kinases can be monitored in lysates of cell lines by real time fluorescent readout using a fluorescent labeled phospho-tyrosine antibody.

## Cell lines

Name	Scirrhous/Non-scirrhous	Histology	Mutation status	response to MTKI		
				Sunitinib	Sorafenib	Pazopanib
4A43	Scirrhous	Signaling cell carcinoma	KRAS G12V			
SB41	Scirrhous	Signaling cell carcinoma	MET amp			
HSC-39	Scirrhous	Signaling cell carcinoma	FGFR2 amp			
HSC-40A	Scirrhous	Signaling cell carcinoma	FGFR2 amp			
HSC-41	Scirrhous	Signaling cell carcinoma	FGFR2 amp			
HSC-40E	Scirrhous	Signaling cell carcinoma	KRAS G12V			
HSC-45	Scirrhous	Signaling cell carcinoma	KRAS G12V			
HSC-38	Scirrhous	Signaling cell carcinoma	MET amp			
HSC-59	Scirrhous	Signaling cell carcinoma	PKCA E545K			
HSC-40	Scirrhous	Signaling cell carcinoma	FGFR2 amp			
KATO-III	Scirrhous	Signaling cell carcinoma	PKCA E542K			
RMS	Non-scirrhous	Moderately-differentiated adenocarcinoma	PKCA E542K			
JH41	Non-scirrhous	Signaling cell carcinoma	PKCA E545K			
NRN1	Non-scirrhous	Adenocarcinoma cell carcinoma	PKCA E545K			
MN28	Non-scirrhous	Moderately-differentiated adenocarcinoma				
MN48	Non-scirrhous	Poorly-differentiated adenocarcinoma	MET amp			
MN7	Non-scirrhous	Well-differentiated adenocarcinoma				
MN74	Non-scirrhous	Moderately-differentiated adenocarcinoma				
NCH07	Non-scirrhous	Well-differentiated adenocarcinoma				
NUC3	Non-scirrhous	Poorty + signal ring component				
NUC4	Non-scirrhous	Poorty + signal ring component				
OCUM1	Non-scirrhous	Signaling cell carcinoma	MET amp			
OKAJMA	Non-scirrhous	Poorly-differentiated adenocarcinoma	MEK1 G55P			
SN-1	Non-scirrhous	Poorly-differentiated adenocarcinoma	KRAS G12D			
SRL18	Non-scirrhous	Poorly-differentiated adenocarcinoma	FGFR2 amp			

Table 1. Results of mutational testing and growth-inhibitory effect of MTKIs (sunitinib, sorafenib, pazopanib).

Sensitive (IC50 values < 1 μmol/L)  
Intermediate (IC50 values > 1 μmol/L)

## Results

Figure 2. Heatmap analysis according to basal tyrosine kinase activity




Figure 3. (A) Volcano plot analysis between sensitive cell lines and non-sensitive cell lines to MTKIs. (B) RT-PCR to detect FGFR2 expression. (C) Copy number analysis on FGFR2




Table 2. Putative kinases relevant to segregate sensitive cell lines

Significant peptides (p<0.01)	Putative kinases contributed to the phosphorylation of peptide
CALM_93_105	ERBB1
IGF1R	IGF1R
INSR	INSR
DYR1A_212_224	DYR1A
ERBB1	ERBB1
EGFR_1165_1177	IGF1R
INSR	INSR
Abl	Abl
EGFR_1190_1202	ERBB1
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
ERBB2	ERBB2
ERBB2	ERBB2
IGF1R	IGF1R
INSR	INSR
ERBB4	ERBB4
IGF1R	IGF1R
INSR	INSR
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
FGFR1	FGFR1
EPHA3	EPHA3
EPHA4	EPHA4
FGFR2	FGFR2
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
IGF1R	IGF1R
INSR	INSR

Significant peptides (p<0.01)	Putative kinases contributed to the phosphorylation of peptide
LAT_194_206	ZAP70
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
IGF1R	IGF1R
INSR	INSR
PDGFRB	PDGFRB
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
IGF1R	IGF1R
INSR	INSR
Src	Src
VEGFR2	VEGFR2

Figure 4. (A) Volcano plot analysis between MET amplified cell lines and MET non-amplified cell lines. (B) Copy number analysis on MET




Table 3. Putative kinases relevant to segregate MET amplified cell lines

Significant peptides (p<0.01)	Putative kinases contributed to the phosphorylation of peptide
EPHA4_689_691	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
ERBB1	ERBB1
ERBB2	ERBB2
ERBB3	ERBB3
ERBB4	ERBB4
EGFR	EGFR
IGF1R	IGF1R
INSR	INSR
ITK	ITK
MET	MET
RET	RET
TEC	TEC
ASK1	ASK1
MAP2K6_MAP2K3	MAP2K6_MAP2K3
MAP2K4_MAP2K7_2	MAP2K4_MAP2K7_2
MEK46	MEK46
rip	rip
VEGFR1	VEGFR1
VEGFR2	VEGFR2
IGF1R	IGF1R
INSR	INSR
ALK	ALK
Pyk	Pyk
PDGFR	PDGFR

Significant peptides (p<0.01)	Putative kinases contributed to the phosphorylation of peptide
RB_804_816	Abl
IGF1R	IGF1R
INSR	INSR
CSF1R	CSF1R
ERBB1	ERBB1
ERBB2	ERBB2
ERBB3	ERBB3
ERBB4	ERBB4
JAK_group	JAK_group
RET	RET
MAP2K6_MAP2K3	MAP2K6_MAP2K3
MAP2K4_MAP2K7_2	MAP2K4_MAP2K7_2
group	group
PDGFRB	PDGFRB
IGF1R	IGF1R
INSR	INSR
VEGFR1	VEGFR1
VEGFR2	VEGFR2
IGF1R	IGF1R
INSR	INSR
ALK	ALK
Pyk	Pyk
PDGFR	PDGFR

Figure 5. Volcano plot analysis between scirrhous cell lines and non-scirrhous cell lines.




Table 4. Putative kinases relevant to segregate scirrhous cell lines

Significant peptides (p<0.01)	Putative kinases contributed to the phosphorylation of peptide
Abl	Abl
ERBB1	ERBB1
ERBB2	ERBB2
ERBB3	ERBB3
ERBB4	ERBB4
IGF1R	IGF1R
INSR	INSR
PDGFRB	PDGFRB
PDGFRB	PDGFRB
Src	Src
ERBB1	ERBB1
IGF1R	IGF1R
INSR	INSR
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
IGF1R	IGF1R
INSR	INSR
ITK	ITK
IGF1R	IGF1R
INSR	INSR
MAP2K6_MAP2K3	MAP2K6_MAP2K3
MAP2K4_MAP2K7_2	MAP2K4_MAP2K7_2
group	group
TEC	TEC
FAK	FAK
IGF1R	IGF1R
INSR	INSR
MEK1	MEK1
INSR	INSR
RET	RET
Src	Src
FAK2	FAK2
IGF1R	IGF1R
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
JAK2	JAK2
EPHA3	EPHA3
EPHA4	EPHA4
EPHA7	EPHA7
IGF1R	IGF1R
INSR	INSR
VEGFR2	VEGFR2
IGF1R	IGF1R
INSR	INSR
SRC	SRC
IGF1R	IGF1R
INSR	INSR

## Conclusions

- Kinase activity profiles reflect the biological complexity of gastric cancer.
- FGFR2 status seems to correlate with sensitivity to sunitinib, pazopanib and sorafenib.
- Tyrosine kinase activity profiling with PamChip® peptide arrays can be utilized to identify molecular targets.
- Pathway analysis needs to be performed to further identify the relevant targets and predictive biomarkers.